

Tabla de referencia VITAPAN 3D-MASTER para VITA In-Ceram® ALUMINA GLASS POWDER

Color VITAPAN 3D-MASTER	1M1	1M2	2L1.5	2L2.5	2M1	2M2	2M3	2R1.5	2R2.5	3L1.5	3L2.5	3M1	3M2	3M3
VITA In-Ceram® ALUMINA GLASS POWDER	AL1	AL1	AL2	AL3	AL1	AL2	AL3	AL2	AL4	AL2	AL3	AL2	AL4	AL3

Color VITAPAN 3D-MASTER	3R1.5	3R2.5	4L1.5	4L2.5	4M1	4M2	4M3	4R1.5	4R2.5	5M1	5M2	5M3
VITA In-Ceram® ALUMINA GLASS POWDER	AL2	AL4	AL2	AL3	AL2	AL4	AL3	AL4	AL4	AL4	AL4	AL3

Tabla de referencia VITAPAN classical para VITA In-Ceram® ALUMINA GLASS POWDER

Color VITAPAN classical	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
VITA In-Ceram® ALUMINA GLASS POWDER	AL1	AL2	AL2	AL3	AL3	AL1	AL2	AL2	AL3	AL3	AL3	AL4	AL4	AL1	AL2	AL3