

CERAMCO II

CERAMCO 3

# TABLA DE HORNEADO PARA HORNOS NEY CENTURION

°F

PROGRAMA EN °F	LO. T °F	DRY. MIN.	RATE °F/MIN.	HI.T °F	HOLD MIN.	COOL MIN.	% VAC	VAC.PULL	VAC.STOP
ULTRA-PAKE	932	3	216	1787	2	2	75%	932	1787
PAINT-O-PAKE PRIMER OPACO	932	3	120	1785	1	2	75%	932	1785
PAINT-O-PAKE SEGUNDO OPACO	932	3	120	1749	1	2	75%	932	1749
EDGE SHOULDER	932	5	126	1769	1	2	75%	932	1769
PRIMER CUERPO	932	5	129	1684	1	2	75%	932	1684
SEGUNDO CUERPO	932	5	129	1675	1	2	75%	932	1675
GLAZE NATURAL	932	3	149	1684	30 SEG.	2	0%	932	1684
ADD-ON	932	3	100	1644	1	2	75%	932	1644
GLAZE DE BAJA TEMPERATURA	932	3	100	1599	1	2	75%	932	1599
MATERIAL REF. MODELO SENCILLO	932	6	50	1850	2	2	0%	932	1850
MATERIAL REF. MODELO MULTIPLE	932	10	50	1850	4	2	0%	932	1850
COLORLOGIC 1° HORNEADO SENCILLO	932	6	75	1680	1	2	75%	932	1680
COLORLOGIC 2° HORNEADO SENCILLO	932	6	100	1670	1	2	75%	932	1670
COLORLOGIC 1° HORNEADO MULTIPLE	932	10	75	1684	1	2	75%	932	1684
COLORLOGIC 2° HORNEADO MULTIPLE	932	10	90	1675	1	2	75%	932	1675
LOW FUSING OPACO Y CUERPO	932	3	100	1590	1	2	75%	932	1590
FINAL TOUCH CUERPO E INCISAL	842	5	99	1301	1	2	75%	842	1301
FINAL TOUCH CARACTERIZADOR Y GLAZE	842	3	99	1274	1.5	2	0%	842	1274

°C

PROGRAMA EN °C	LO. T °C	DRY. MIN.	RATE °C/MIN.	HI.T °C	HOLD MIN.	COOL MIN.	% VAC	VAC.PULL	VAC.STOP
ULTRA-PAKE	500	3	120	975	2	2	75%	500	975
PAINT-O-PAKE PRIMER OPACO	500	3	67	974	1	2	75%	500	974
PAINT-O-PAKE SEGUNDO OPACO	500	3	67	954	1	2	75%	500	954
EDGE SHOULDER	500	5	70	965	1	2	75%	500	965
PRIMER CUERPO	500	5	72	918	1	2	75%	500	918
SEGUNDO CUERPO	500	5	72	913	1	2	75%	500	913
GLAZE NATURAL	500	3	83	918	30 SEG.	2	0%	500	918
ADD-ON	500	3	56	896	1	2	75%	500	896
GLAZE DE BAJA TEMPERATURA	500	3	56	871	1	2	75%	500	871
MATERIAL REF. MODELO SENCILLO	500	6	28	1010	2	2	0%	500	1010
MATERIAL REF. MODELO MULTIPLE	500	10	28	1010	4	2	0%	500	1010
COLORLOGIC 1° HORNEADO SENCILLO	500	6	42	916	1	2	75%	500	916
COLORLOGIC 2° HORNEADO SENCILLO	500	6	56	910	1	2	75%	500	910
COLORLOGIC 1° HORNEADO MULTIPLE	500	10	42	918	1	2	75%	500	918
COLORLOGIC 2° HORNEADO MULTIPLE	500	10	50	913	1	2	75%	500	913
LOW FUSING OPACO Y CUERPO	500	3	56	866	1	2	75%	500	866
FINAL TOUCH CUERPO E INCISAL	450	5	55	705	1	2	75%	450	705
FINAL TOUCH CARACTERIZADOR Y GLAZE	450	3	55	690	1.5	2	0%	450	690

# MEDIDENTAL

PALENQUE # 456 COL. NARVARTE

MEXICO D. F. C. P. 03600

TEL: 55-59-48-65 FAX: 55-75-43-03

E-mail: federicojassan@hotmail.com